



GCM2 Antibody, Biotin conjugated

Product Code	CSB-PA009324LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75603
Immunogen	Recombinant Human Chorion-specific transcription factor GCMb protein (156-451AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chorion-specific transcription factor GCMb (hGCMb) (GCM motif protein 2) (Glial cells missing homolog 2), GCM2, GCMB
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	GCM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.