




GDF11 Antibody

Product Code	CSB-PA009344LA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BG54
Immunogen	Recombinant Pig Growth differentiation factor 11 protein (296-404AA)
Raised In	Rabbit
Species Reactivity	Pig, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDF11
Immunogen Species	Pig
Research Area	Others
Target Names	GDF11

Image



Western Blot

Positive WB detected in: Mouse kidney tissue,
Mouse brain tissue

All lanes: GDF11 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 46 kDa

Observed band size: 46 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.