



GDF11 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA009344LC01PI |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BG54 |
| Immunogen | Recombinant Pig Growth differentiation factor 11 protein (296-404AA) |
| Raised In | Rabbit |
| Species Reactivity | Pig |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | GDF11 |
| Immunogen Species | Pig |
| Research Area | Others |
| Target Names | GDF11 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |