



Gamma-glutamyl hydrolase Antibody, Biotin conjugated

Product Code	CSB-PA009389LD01GGV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P93164
Immunogen	Recombinant Glycine max Gamma-glutamyl hydrolase protein (22-342AA)
Raised In	Rabbit
Species Reactivity	Glycine max
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-glutamyl hydrolase (EC 3.4.19.9) (Conjugase) (GH) (Gamma-Glu-X carboxypeptidase)
Immunogen Species	Glycine max
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.