



GJA8 Antibody, Biotin conjugated

Product Code	CSB-PA009449LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48165
Immunogen	Recombinant Human Gap junction alpha-8 protein (286-432AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gap junction alpha-8 protein (Connexin-50) (Cx50) (Lens fiber protein MP70), GJA8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GJA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.