





GNAQ Antibody

Product Code	CSB-PA009594LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50148
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(q) subunit alpha protein (1-359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein G(q) subunit alpha (Guanine nucleotide-binding protein alpha-q), GNAQ, GAQ
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	GNAQ

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA009594LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA009594LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.