



GPR85 Antibody, HRP conjugated

Product Code	CSB-PA009837LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60893
Immunogen	Recombinant Human Probable G-protein coupled receptor 85 protein (229-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 85 (Super conserved receptor expressed in brain 2), GPR85, SREB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GPR85
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.