



GPR85 Antibody, HRP conjugated

Product Code	CSB-PA009837LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60893
Immunogen	Recombinant Human Probable G-protein coupled receptor 85 protein (229-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 85 (Super conserved receptor expressed in brain 2), GPR85, SREB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GPR85