



HMBS Antibody, Biotin conjugated

Product Code	CSB-PA010524LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08397
Immunogen	Recombinant Human Porphobilinogen deaminase protein (283-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Porphobilinogen deaminase (PBG-D) (EC 2.5.1.61) (Hydroxymethylbilane synthase) (HMBS) (Pre-uroporphyrinogen synthase), HMBS, PBGD UPS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HMBS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.