




HOXA5 Antibody

Product Code	CSB-PA010655LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20719
Immunogen	Recombinant Human Homeobox protein Hox-A5 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-A5 (Homeobox protein Hox-1C), HOXA5, HOX1C
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXA5

Image



Western Blot

Positive WB detected in: SH-SY5Y whole cell lysate, HepG2 whole cell lysate, A549 whole cell lysate, Rat heart tissue


All lanes: HOXA5 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 49 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA010655LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.