



# HSD3B1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA010781LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14060
<b>Immunogen</b>	Recombinant Human 3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 1 protein (127-267AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 1 (3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type I) (3-beta-HSD I) (Trophoblast antigen FDO161G) [Includes: 3-beta-hydroxy-Delta(5)-steroid dehydrogenase (EC 1.1.1.145) (3-beta-hydroxy-5-ene steroid dehydrogenase) (Progesterone reductase); Steroid Delta-isomerase (EC 5.3.3.1) (Delta-5-3-ketosteroid isomerase)], HSD3B1, 3BH HSDB3A
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	HSD3B1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.