



HSPA4 Antibody, FITC conjugated

Product Code	CSB-PA010825LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34932
Immunogen	Recombinant Human Heat shock 70 kDa protein 4 protein (1-148AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock 70 kDa protein 4 (HSP70RY) (Heat shock 70-related protein APG-2), HSPA4, APG2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPA4