




HTR3A Antibody

Product Code	CSB-PA010890LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46098
Immunogen	Recombinant Human 5-hydroxytryptamine receptor 3A protein (322-455AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-hydroxytryptamine receptor 3A (5-HT3-A) (5-HT3A) (5-hydroxytryptamine receptor 3) (5-HT-3) (5-HT3R) (Serotonin receptor 3A) (Serotonin-gated ion channel receptor), HTR3A, 5HT3R HTR3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HTR3A

Image



Western Blot

Positive WB detected in: Rat lung tissue, Mouse lung tissue


All lanes: HTR3A antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56, 59, 54, 60 kDa

Observed band size: 56 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA010890LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.