



# HTR3A Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA010890LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P46098
<b>Immunogen</b>	Recombinant Human 5-hydroxytryptamine receptor 3A protein (322-455AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	5-hydroxytryptamine receptor 3A (5-HT3-A) (5-HT3A) (5-hydroxytryptamine receptor 3) (5-HT-3) (5-HT3R) (Serotonin receptor 3A) (Serotonin-gated ion channel receptor), HTR3A, 5HT3R HTR3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	HTR3A
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.