



HTRA2 Antibody, FITC conjugated

Product Code	CSB-PA010902EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43464
Immunogen	Recombinant Human Serine protease HTRA2, mitochondrial protein (334-458AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease HTRA2, mitochondrial (EC 3.4.21.108) (High temperature requirement protein A2) (HtrA2) (Omi stress-regulated endoprotease) (Serine protease 25) (Serine proteinase OMI), HTRA2, OMI PRSS25
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HTRA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.