




IFNG Antibody

Product Code	CSB-PA011050YA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07353
Immunogen	Recombinant Bovine Interferon gamma protein (24-166AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma (BoIFNG) (IFN-gamma), IFNG
Immunogen Species	Bos taurus (Bovine)
Research Area	Immunology
Target Names	IFNG

Image



Western Blot
 Positive WB detected in Recombinant protein
 All lanes: IFNG antibody at 3µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 predicted band size: 20 kDa
 observed band size: 20 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.