



IFNG Antibody, Biotin conjugated

Product Code	CSB-PA011050YD01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07353
Immunogen	Recombinant Bovine Interferon gamma protein (24-166AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon gamma (BoIFNG) (IFN-gamma), IFNG
Immunogen Species	Bos taurus (Bovine)
Research Area	Immunology
Target Names	IFNG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.