



IL12RB1 Antibody, Biotin conjugated

Product Code	CSB-PA011588LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42701
Immunogen	Recombinant Human Interleukin-12 receptor subunit beta-1 protein (119-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 receptor subunit beta-1 (IL-12 receptor subunit beta-1) (IL-12R subunit beta-1) (IL-12R-beta-1) (IL-12RB1) (IL-12 receptor beta component) (CD antigen CD212), IL12RB1, IL12R IL12RB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12RB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.