



IL12RB1 Antibody, Biotin conjugated

Product Code	CSB-PA011588LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42701
Immunogen	Recombinant Human Interleukin-12 receptor subunit beta-1 protein (119-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 receptor subunit beta-1 (IL-12 receptor subunit beta-1) (IL-12R subunit beta-1) (IL-12R-beta-1) (IL-12RB1) (IL-12 receptor beta component) (CD antigen CD212), IL12RB1, IL12R IL12RB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12RB1