




IRF6 Antibody

Product Code	CSB-PA011821LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14896
Immunogen	Recombinant Human Interferon regulatory factor 6 protein (43-285AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon regulatory factor 6 (IRF-6), IRF6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRF6

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate, A549 whole cell lysate, MCF-7 whole cell lysate, , Mouse liver tissue


All lanes: IRF6 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 54, 42 kDa

Observed band size: 54 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA011821LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.