



# IRF6 Antibody, HRP conjugated

Product Code	CSB-PA011821LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14896
Immunogen	Recombinant Human Interferon regulatory factor 6 protein (43-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon regulatory factor 6 (IRF-6), IRF6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRF6