



ITGB2 Antibody, Biotin conjugated

Product Code	CSB-PA011884LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05107
Immunogen	Recombinant Human Integrin beta-2 protein (37-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin beta-2 (Cell surface adhesion glycoproteins LFA-1/CR3/p150,95 subunit beta) (Complement receptor C3 subunit beta) (CD antigen CD18), ITGB2, CD18 MFI7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.