



KCNA1 Antibody, HRP conjugated

Product Code	CSB-PA012005LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q09470
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily A member 1 protein (7-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily A member 1 (Voltage-gated K(+) channel HuK1) (Voltage-gated potassium channel HBK1) (Voltage-gated potassium channel subunit Kv1.1), KCNA1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.