



KCNC1 Antibody, Biotin conjugated

Product Code	CSB-PA012018LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48547
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily C member 1 protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily C member 1 (NGK2) (Voltage-gated potassium channel subunit Kv3.1) (Voltage-gated potassium channel subunit Kv4), KCNC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.