




KCNQ4 Antibody

Product Code	CSB-PA012092LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56696
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily KQT member 4 protein (401-516AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily KQT member 4 (KQT-like 4) (Potassium channel subunit alpha KvLQT4) (Voltage-gated potassium channel subunit Kv7.4), KCNQ4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNQ4

Image



Immunofluorescent analysis of U251 cells using CSB-PA012092LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.