



KLKB1 Antibody, HRP conjugated

Product Code	CSB-PA012461LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03952
Immunogen	Recombinant Human Plasma kallikrein protein (503-612AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasma kallikrein (EC 3.4.21.34) (Fletcher factor) (Kininogenin) (Plasma prekallikrein) (PKK) [Cleaved into: Plasma kallikrein heavy chain; Plasma kallikrein light chain], KLKB1, KLK3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	KLKB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.