



KRT7 Antibody, Biotin conjugated

Product Code	CSB-PA012564LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08729
Immunogen	Recombinant Human Keratin, type II cytoskeletal 7 protein (2-469AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 7 (Cytokeratin-7) (CK-7) (Keratin-7) (K7) (Sarcolelectin) (Type-II keratin Kb7), KRT7, SCL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.