



LDHC Antibody, FITC conjugated

Product Code	CSB-PA012844LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07864
Immunogen	Recombinant Human L-lactate dehydrogenase C chain protein (2-332AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-lactate dehydrogenase C chain (LDH-C) (EC 1.1.1.27) (Cancer/testis antigen 32) (CT32) (LDH testis subunit) (LDH-X), LDHC, LDH3 LDHX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LDHC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.