




LEFTY1 Antibody

Product Code	CSB-PA012857LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75610
Immunogen	Recombinant Human Left-right determination factor 1 protein (223-260AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Left-right determination factor 1 (Left-right determination factor B) (Protein lefty-1) (Protein lefty-B), LEFTY1, LEFTB LEFTYB
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	LEFTY1

Image



Western Blot

Positive WB detected in: Rat spleen tissue, Mouse heart tissue, Mouse liver tissue, Mouse brain tissue

All lanes: LEFTY1 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41 kDa

Observed band size: 41 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.