



LEFTY1 Antibody, Biotin conjugated

Product Code	CSB-PA012857LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75610
Immunogen	Recombinant Human Left-right determination factor 1 protein (223-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Left-right determination factor 1 (Left-right determination factor B) (Protein lefty-1) (Protein lefty-B), LEFTY1, LEFTB LEFTYB
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	LEFTY1