



MAEA Antibody, Biotin conjugated

Product Code	CSB-PA013310LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L5Y9
Immunogen	Recombinant Human Macrophage erythroblast attacher protein (299-356AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage erythroblast attacher (Cell proliferation-inducing gene 5 protein) (Erythroblast macrophage protein) (Human lung cancer oncogene 10 protein) (HLC-10), MAEA, EMP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MAEA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.