



MAP2K1 Antibody, HRP conjugated

Product Code	CSB-PA013409LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02750
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 1 protein (2-79AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity mitogen-activated protein kinase kinase 1 (MAP kinase kinase 1) (MAPKK 1) (MKK1) (EC 2.7.12.2) (ERK activator kinase 1) (MAPK/ERK kinase 1) (MEK 1), MAP2K1, MEK1 PRKMK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP2K1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.