





MAPK11 Antibody

Product Code	CSB-PA013450LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15759
Immunogen	Recombinant Human Mitogen-activated protein kinase 11 protein (310-362AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 11 (MAP kinase 11) (MAPK 11) (EC 2.7.11.24) (Mitogen-activated protein kinase p38 beta) (MAP kinase p38 beta) (p38b) (Stress-activated protein kinase 2b) (SAPK2b) (p38-2), MAPK11, PRKM11 SAPK2 SAPK2B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAPK11

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA013450LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA013450LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.