



MAPK12 Antibody, HRP conjugated

Product Code	CSB-PA013451LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53778
Immunogen	Recombinant Human Mitogen-activated protein kinase 12 protein (314-367AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 12 (MAP kinase 12) (MAPK 12) (EC 2.7.11.24) (Extracellular signal-regulated kinase 6) (ERK-6) (Mitogen-activated protein kinase p38 gamma) (MAP kinase p38 gamma) (Stress-activated protein kinase 3), MAPK12, ERK6 SAPK3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAPK12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.