



MLLT4 Antibody, Biotin conjugated

Product Code	CSB-PA014636LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55196
Immunogen	Recombinant Human Afadin protein (1663-1809AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Afadin (ALL1-fused gene from chromosome 6 protein) (Protein AF-6) (Afadin adherens junction formation factor), AFDN, AF6 MLLT4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	AFDN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.