



MLPH Antibody, HRP conjugated

Product Code	CSB-PA014640LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV36
Immunogen	Recombinant Human Melanophilin protein (310-539AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanophilin (Exophilin-3) (Slp homolog lacking C2 domains a) (Slac2-a) (Synaptotagmin-like protein 2a), MLPH, SLAC2A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MLPH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.