



MSC Antibody, HRP conjugated

Product Code	CSB-PA015025LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60682
Immunogen	Recombinant Human Musculin protein (3-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Musculin (Activated B-cell factor 1) (ABF-1) (Class A basic helix-loop-helix protein 22) (bHLHa22), MSC, ABF1 BHLHA22
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MSC