



# MTRF1 Antibody, HRP conjugated

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA015197LB01HU  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                       |
| Uniprot No.         | O75570  |
| Immunogen           | Recombinant Human Peptide chain release factor 1, mitochondrial protein (285-433AA) |
| Raised In           | Rabbit  |
| Species Reactivity  | Human   |
| Tested Applications | ELISA   |
| Form                | Liquid  |
| Conjugate           | HRP   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4    |
| Purification Method | >95%, Protein G purified  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Peptide chain release factor 1, mitochondrial (MRF-1) (MtRF-1), MTRF1               |
| Immunogen Species   | Homo sapiens (Human)  |
| Research Area       | Epigenetics and Nuclear Signaling   |
| Target Names        | MTRF1   |