



MTRF1 Antibody, HRP conjugated

Product Code	CSB-PA015197LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75570
Immunogen	Recombinant Human Peptide chain release factor 1, mitochondrial protein (285-433AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptide chain release factor 1, mitochondrial (MRF-1) (MtRF-1), MTRF1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MTRF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.