



NCF1 Antibody, HRP conjugated

Product Code	CSB-PA015525LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14598
Immunogen	Recombinant Human Neutrophil cytosol factor 1 protein (63-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutrophil cytosol factor 1 (NCF-1) (47 kDa autosomal chronic granulomatous disease protein) (47 kDa neutrophil oxidase factor) (NCF-47K) (Neutrophil NADPH oxidase factor 1) (Nox organizer 2) (Nox-organizing protein 2) (SH3 and PX domain-containing protein 1A) (p47-phox), NCF1, NOXO2 SH3PXD1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NCF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.