



NLRX1 Antibody, Biotin conjugated

Product Code	CSB-PA015878LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UT6
Immunogen	Recombinant Human NLR family member X1 protein (45-205AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NLR family member X1 (Caterpillar protein 11.3) (CLR11.3) (Nucleotide-binding oligomerization domain protein 26) (Nucleotide-binding oligomerization domain protein 5) (Nucleotide-binding oligomerization domain protein 9), NLRX1, NOD26 NOD5 NOD9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.