



NR5A2 Antibody, Biotin conjugated

Product Code	CSB-PA016066LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00482
Immunogen	Recombinant Human Nuclear receptor subfamily 5 group A member 2 protein (182-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 5 group A member 2 (Alpha-1-fetoprotein transcription factor) (B1-binding factor) (hB1F) (CYP7A promoter-binding factor) (Hepatocytic transcription factor) (Liver receptor homolog 1) (LRH-1), NR5A2, B1F CPF FTF
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	NR5A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.