



NRAS Antibody, HRP conjugated

Product Code	CSB-PA016070LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01111
Immunogen	Recombinant Human GTPase NRas protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase NRas (Transforming protein N-Ras), NRAS, HRAS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NRAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.