



NUAK1 Antibody

Product Code	CSB-PA016139LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60285
Immunogen	Recombinant Human NUAK family SNF1-like kinase 1 protein (442-573AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NUAK family SNF1-like kinase 1 (EC 2.7.11.1) (AMPK-related protein kinase 5) (ARK5) (Omphalocele kinase 1), NUAK1, ARK5 KIAA0537 OMPHK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUAK1