



NUAK1 Antibody, FITC conjugated

Product Code	CSB-PA016139LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60285
Immunogen	Recombinant Human NUAK family SNF1-like kinase 1 protein (442-573AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NUAK family SNF1-like kinase 1 (EC 2.7.11.1) (AMPK-related protein kinase 5) (ARK5) (Omphalocele kinase 1), NUAK1, ARK5 KIAA0537 OMPHK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUAK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.