



NUDT1 Antibody, Biotin conjugated

Product Code	CSB-PA016154LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36639
Immunogen	Recombinant Human 7,8-dihydro-8-oxoguanine triphosphatase protein (49-82AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	7,8-dihydro-8-oxoguanine triphosphatase (EC 3.6.1.55) (2-hydroxy-dATP diphosphatase) (EC 3.6.1.56) (8-oxo-dGTPase) (Nucleoside diphosphate-linked moiety X motif 1) (Nudix motif 1), NUDT1, MTH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NUDT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.