



NUP107 Antibody, FITC conjugated

Product Code	CSB-PA016188LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57740
Immunogen	Recombinant Human Nuclear pore complex protein Nup107 protein (1-66AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup107 (107 kDa nucleoporin) (Nucleoporin Nup107), NUP107
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NUP107
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.