




Ocm Antibody

Product Code	CSB-PA016264HA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51879
Immunogen	Recombinant Mouse Oncomodulin protein (2-109AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oncomodulin (OM) (Parvalbumin beta), Ocm
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Ocm

Image



Western Blot

Positive WB detected in: Mouse heart tissue

All lanes: Ocm antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 13 kDa

Observed band size: 13 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.