



Ocm Antibody, FITC conjugated

Product Code	CSB-PA016264HC01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51879
Immunogen	Recombinant Mouse Oncomodulin protein (2-109AA)
Raised In	Rabbit
Species Reactivity	Mouse
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oncomodulin (OM) (Parvalbumin beta), Ocm
Immunogen Species	Mus musculus (Mouse)
Research Area	Neuroscience
Target Names	Ocm
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.