



OXT Antibody, Biotin conjugated

Product Code	CSB-PA017315LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01178
Immunogen	Recombinant Human Oxytocin-neurophysin 1 protein (27-123AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxytocin-neurophysin 1 (OT-NPI) [Cleaved into: Oxytocin (Ocytocin); Neurophysin 1], OXT, OT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OXT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.